



Pediatric Examination Table: Model 5900

CHAIR SPECIFICATIONS

- Overall height: 44.5" (113cm)
- Work surface height: 35" (88.9cm)
- Width: 48" (121.9cm)
- Depth: 25" (63.5cm)
- Drawers: 16.5" W x 17" D x 4" H (42cm x 43cm x 10cm)
- Overall door opening: 40.25" W x 17.25" H (102cm x 44cm)
- Weight capacity: 50lb (22.7kg)

KEY FEATURES

- Ultra comfort seamless top with antimicrobial guard
- Heavy duty, all-metal storage drawers
- Heavy duty roller bearing drawer glides
- All metal cabinet doors with commercial grade hinges
- Padded front rail for infant safety
- 18" (45.7cm) paper roll holder
- Digital scale operated by 6 D-cell batteries (not included) with auto shut-off feature
- SECA® stainless steel measuring rule with paddles
- Weighs in both pounds/ounces and kilograms
- Tare weight capability
- 3 year limited capability

ACCESSORIES

- 091 – KG scale only
- 092 – Pounds scale only
- 365 – I.V. pole for 5200 series tables (*factory installed only*)
- 502 – Double door lock (*factory installed only*)
- 505 – Single drawer lock (*factory installed only*)

COLOR:

